



## UltraClean™ Agarose, Molecular Biology Grade

*For the electrophoretic separation of nucleic acids*

**DESCRIPTION** UltraClean™ Agarose, Molecular Biology Grade, is a high quality agarose recommended for the routine separation and recovery of nucleic acids up from 0.2 to 50 kb. Each lot of UltraClean™ Agarose, Molecular Biology Grade, is tested and certified to be free of DNase and RNase activity.



## UltraClean™ MS-8 Agarose

*For high resolution of small DNA fragments and PCR products*

**DESCRIPTION** UltraClean™ MS-8 Agarose is for molecular screening that improves resolution of small DNA fragments and PCR products from 50 – 1000 bp. This performance agarose features improved clarity, higher gel strength and is DNase and RNase free.



## UltraClean™ Forensic Agarose

*For the electrophoretic separation of nucleic acids >1000 bp*

**DESCRIPTION** UltraClean™ Forensic Agarose is ideal for laboratories performing forensic testing, paternity determination, cell line verification and tissue typing. This agarose is certified to be DNA-free, DNase-free and RNase-free and meets all requirements for DNA identity applications.



## UltraClean™ Low Melt Agarose

*A low-melting agarose for in-gel reactions*

**DESCRIPTION** UltraClean™ Low Melt Agarose is ideal for in-gel DNA ligation reactions. This agarose melts at 65°C and can resolve a broad range of DNA fragments from 0.5 to 20 kb.



## UltraClean™ Low Melt Sieve Agarose

*For the electrophoretic separation of nucleic acids 10 – 1000 bp*

**DESCRIPTION** UltraClean™ Low Melt Sieve Agarose is a premium agarose useful for separation and recovery of low molecular weight DNA fragments of 10 to 1000 bp from agarose. An economical replacement for NuSieve® agarose.

UltraClean™ Agarose, Molecular Biology Grade	
Catalog No.	Quantity
15003-50	50 g
15003-100	100 g
15003-500	500 g
15003-1000	1 kg
UltraClean™ MS-8 Agarose	
Catalog No.	Quantity
15515-50	50 g
15515-100	100 g
15515-500	500 g

UltraClean™ Forensic Agarose	
Catalog No.	Quantity
15505-50	50 g
15505-100	100 g
15505-500	500 g
UltraClean™ Low Melt Agarose	
Catalog No.	Quantity
15005-50	50 g
15005-100	100 g
15005-500	500 g
UltraClean™ Low Melt Sieve Agarose	
Catalog No.	Quantity
15004-50	50 g
15004-100	100 g
15004-500	500 g