

BOECO MAGNETIC STIRRER WITH HOTPLATE MSH 420

Motor drive Magnetic stirrer with all in one glass-ceramic hot plate for simultaneously stirring and heating of chemical reagents.

The Glass-ceramic top plate is easy to clean, chemical-resistant, corrosion-resistant and scratch resistant.

The MSH 420 provides LED Indicators for Power, Heating and Stirring.

Specification

| | |
|-------------------------|---|
| Speed range: | 60 - 1110 1/min |
| Temperature regulation: | ambient temperature to 450 C°. |
| Max. stirring volume: | up to 15,0 L. |
| Hot Platform size: | 200 x 200 mm |
| Hot Platform material: | glass-ceramic |
| "Hot Top" indicator: | a "Hot" light illuminates if the temperature on the hot plate exceeds 55°C |
| Over-Heat protection: | The heating power will cutoff in abnormal situations by limiting top plate temperature to 550 °C, and light flash |
| Overall size: | 320 (W) x 210 (D) x 105 (H) mm |
| Heating power: | 750 Watts |
| Weight: | 3,6 kg |



MSH 420

| Code | Description |
|-------------|--|
| BOE 8054000 | Magnetic stirrer with hotplate MSH 420, 220V |
| BOE 8054100 | Magnetic stirrer with hotplate MSH 420, 110V |

BOECO SUNFLOWER MINI-SHAKER 3D

Sunflower mini Shaker 3D with adjustable speed and constant angle of mixing. Very small and ideal for private diagnostics

Incl. black non-slip rubber mat

Specification

| | |
|-----------------|---------------------------------------|
| Kind of motion: | three-dimensional orbital movement |
| Angle: | 7 ° |
| Speed range: | 10 - 30 rpm |
| Speed display: | scale (min/max) |
| Platform size: | 200 x 200 mm |
| Overall size: | 85 x 150 x 80 mm |
| Max Loading: | 500 g |
| Weight: | 1,0 kg |
| Power: | 12 V. 300 mA (External power adapter) |



3D

| Code | Description |
|-------------|---|
| BOE 8058000 | Sunflower Mini Shaker 3D, AC Adapter 220V |
| BOE 8058100 | Sunflower Mini Shaker 3D, AC Adapter 110V |
| BOE 8050100 | TP-26 tube attachment, for tubes 13 x 75 mm |

BOECO ROTATOR RS-24

Mini-rotator with vertical rotation of the platform for soft mixing of test tubes, microtubes and vacutainers.

Incl. Standard platform PRS-22 for 22 tubes up to 15 mm diam.

Specification

| | |
|----------------|---------------------------------------|
| Speed range: | 5 - 30 1/min. |
| Digital Timer: | from 1 min to 23 H 59 min |
| Overall size: | 310 x 40 x 175 mm |
| Weight: | 1,5 kg |
| Power: | 12 V. 500 mA (External power adapter) |

| Code | Description |
|-------------|---|
| BOE 8024000 | Rotator RS-24, with platform PRS-22 AC adapter 220V 50/60 Hz |
| BOE 8024100 | Rotator RS-24, with platform PRS-22 AC adapter 110V 60 Hz |
| BOE 8073001 | Platform PRS-4/12 (4x30 mm diam and 12x15 mm) |



RS-24